






 **TAEWOONG**
NITI-S™




Gastrointestinal Self-expandable Metal Stent

Duodenal & Colonic Stent

Biliary Stent


D



for malignant obstruction

- **Unfixed cell with weaving construction**
 - low foreshortening for accurate positioning
 - optimal combination of radial and axial force to maintain luminal patency in tortuous anatomy

S



"Fully and partially covered available"

for benign and malignant obstruction

- **Silicone coated design with a large trunk**
 - to prevent the risk of tumor ingrowth
- **Retrieval string facilitates safe and smooth removal**


Esophageal Stent

COMVI™



for malignant obstruction

- **Triple layered construction**
 - PTFE membrane prevents the risk of tissue invasion
 - outer wire mesh prevents the risk of migration
 - unfixed cell structure enables the stent to fit in tortuous anatomy



COMVI™ (Flare)



for malignant obstruction

- **Triple layered construction**
 - uncovered proximal flare end to reduce migration
 - large cell structured body with thick wire to reduce the risk of fracture and enhance radial force
- **Large diameter (up to 26mm)** is loaded into a 10.5Fr delivery system

Duodenal & Colonic Stent

Ordering Information

Stent Type	OTW			TTS			Delivery Length (cm)	
	Code	Diameter (mm)	Length (cm)	Code	Diameter (mm)	Length (cm)	OTW	TTS
Pyloric / Duodenal Stent								
S (Covered)	PS**##F	18, 20, 22, 24, 26, 28	6, 8, 10, 12, 14, 15	PST**##F	18, 20, 22	6, 8, 10, 12, 14, 15	135	180
S (Both bare)	PS**##B			PST**##B				
S (Proximal bare)	PS**##H			PST**##H				
D	PD**##			PDT**##	18, 20, 22, 24			
COMVI™	PC**##BA		6, 8, 10, 12	PCT**##BA	18, 20, 22	6, 8, 10, 12		
COMVI™ (Flare)				PCT**##P	18, 20, 22, 24, 26	6, 8, 10, 12		180
Enteral Colonic Stent								
S (Covered)	CS**##F	18, 20, 22, 24, 26, 28	6, 8, 10, 12, 14, 15	CST**##F	18, 20, 22	6, 8, 10, 12, 14, 15	70	220
S (Both bare)	CS**##B			CST**##B				
S (Distal bare)	CS**##H			CST**##H				
D	CD**##	18,20,22,24,26,28,30		CDT**##	18,20,22,24,26,28			
COMVI™	CC**##BA	18,20,22,24,26,28	6, 8, 10, 12	CCT**##BA	18, 20, 22	6, 8, 10, 12		
COMVI™ (Flare)				CCT**##P	18, 20, 22, 24, 26	6, 8, 10, 12		220